

# WTSoftwareUtilitiesShare

## Wind Tunnel Analysis Software Package

### Repository Structure and Usage Guide

Analysis Documentation

December 10, 2025

## Contents

<b>1</b>	<b>Repository Overview</b>	<b>2</b>
1.1	Repository Structure . . . . .	2
<b>2</b>	<b>Main Workflows and Example Files</b>	<b>3</b>
2.1	Example_PointConc_Analysis_200325_stepByStep (Learning) . . . . .	3
2.1.1	Processing Steps . . . . .	3
2.2	Example_Flow_Analysis (Flow Field Analysis) . . . . .	5
2.3	Example_PointConc_Analysis_200325_Analysis_deploy . . . . .	5
<b>3</b>	<b>Physical Equations for Concentration Analysis</b>	<b>6</b>
3.1	Model Scale Mass Flow Rate . . . . .	6
3.2	Net Concentration . . . . .	6
3.3	Dimensionless Concentration ( $C^*$ ) . . . . .	6
3.4	Full-Scale Concentration . . . . .	7
<b>4</b>	<b>Parameter File Structure</b>	<b>9</b>
4.1	Example CSV Format . . . . .	9
<b>5</b>	<b>Software Usage</b>	<b>10</b>
5.1	Quick Start for Students (Jupyter Server) . . . . .	10
5.2	Local Installation (Windows) . . . . .	10
5.3	Manual Installation (Linux/Mac) . . . . .	10
<b>6</b>	<b>Conclusion</b>	<b>10</b>

# 1 Repository Overview

The `WTSoftwareUtilitiesShare` repository provides a comprehensive Python package for wind tunnel measurement analysis, focusing on point concentration measurements and flow analysis. The software enables automated data processing, dimensionalization, and visualization of wind tunnel experimental data.

## 1.1 Repository Structure

- **Main Package:** `windtunnel/` – Core Python package containing analysis modules
  - `concentration/` – Point and puff concentration analysis
  - `flow/` – Flow statistics and turbulence analysis
  - `plots/` – Visualization tools
  - `grid/` – Grid manipulation tools
- **Jupyter Notebooks:**
  - `Example_PointConc_Analysis_200325_stepByStep.ipynb` – Tutorial notebook
  - `Example_PointConc_Analysis_200325_Analysis_deploy.ipynb` – Deployment version
  - `Example_PointConc_Analysis_200325_Analysis.ipynb` – Full analysis workflow
  - `Example_Flow_Analysis.ipynb` – Flow field analysis
- **Data Structure:**
  - `Data/InputData/` – Raw measurement time series
  - `Data/ParameterFiles/` – CSV files with ambient conditions
  - `Results/` – Processed output data
  - `ExampleData/` – Example datasets for testing
- **Deployment:**
  - `deploy_config.py` – Setup functions for Jupyter server deployment
  - `windows_deploy/` – Windows batch scripts for environment setup
  - `requirements.txt` – Python dependencies

## 2 Main Workflows and Example Files

### 2.1 Example\_PointConc\_Analysis\_200325\_stepByStep (Learning)

This notebook provides a detailed step-by-step walkthrough of the concentration analysis workflow. It is designed for learning and understanding the underlying processes.

#### 2.1.1 Processing Steps

##### 1. Data Import and Initialization

- Import `windtunnel` package and required libraries
- Define paths to input data and output directories
- Specify measurement file prefixes (namelist)
- Define CSV parameter file location

##### 2. Read Ambient Conditions

- Load ambient conditions from CSV file for each time series
- Extract source location ( $x_{\text{source}}, y_{\text{source}}, z_{\text{source}}$ )
- Extract measurement location ( $x_{\text{measure}}, y_{\text{measure}}, z_{\text{measure}}$ )
- Read environmental parameters (pressure, temperature)
- Read mass flow controller settings and calibration data
- Read scaling information (model scale, reference lengths)
- Read tracer gas properties (molecular weight, gas factor)
- Read full-scale transformation parameters

##### 3. Initialize PointConcentration Objects

- Create dictionary structure for all measurement files
- Load raw time series data using `PointConcentration.from_file()`
- Each file contains: time, wtref, slow\_FID, fast\_FID, open\_rate
- Store data in nested dictionaries: `conc_ts[name][file]`

##### 4. Set Ambient Conditions

- Apply ambient conditions to each `PointConcentration` object
- Store source and measurement positions
- Store pressure ( $p$ ) and temperature ( $T$ )
- Store mass flow controller settings
- Calculate relative measurement position:  $x = x_{\text{measure}} - x_{\text{source}}$

##### 5. Set Scaling Information

- Define model to reality scale factor
- Set reference length  $L_{\text{ref}}$  (model scale)
- Set reference height  $H_{\text{ref}}$
- Define scaling factor for wtref adjustment
- Calculate full-scale reference length:  $L_{\text{ref,fs}} = \text{scale} \cdot L_{\text{ref}}$

## 6. Set Tracer Gas Information

- Define gas name (e.g., C12, N<sub>2</sub>)
- Set molecular weight  $M$  [g/mol]
- Set gas correction factor (for different gases)

## 7. Set Full-Scale Information

- Define full-scale reference wind speed  $U_{\text{ref,fs}}$  [m/s]
- Define full-scale mass flow rate  $\dot{Q}_{\text{fs}}$  [kg/s]
- Set full-scale ambient temperature  $T_{\text{fs}}$  [°C]
- Set full-scale ambient pressure  $p_{\text{fs}}$  [Pa]

## 8. Temperature Conversion

- Convert ambient temperature to Kelvin:  $T_K = T + 273.15$
- Convert full-scale temperature:  $T_{\text{fs,K}} = T_{\text{fs}} + 273.15$
- Store standard conditions:  $T_{\text{std,K}} = 273.15$  K,  $p_{\text{std}} = 101325$  Pa

## 9. Calculate Mean Wind Tunnel Reference Speed

- Calculate mean of wtref time series
- Apply scaling factor:  $\bar{U}_{\text{ref,ms}} = \text{scaling\_factor} \cdot \text{mean}(\text{wtref})$

## 10. Calculate Model Scale Mass Flow Rate

- Account for controller settings and open\_rate
- Apply gas correction factor
- Correct for ambient conditions (temperature and pressure)
- Apply calibration (if specified)
- **See Equation 1 for details**

## 11. Calculate Net Concentration

- Subtract background concentration from total concentration
- $C_{\text{net}} = \text{fast\_FID} - \text{slow\_FID}$  [ppmV]

## 12. Calculate Dimensionless Concentration ( $C^*$ )

- Normalize concentration by flow rate and reference wind speed
- **See Equation 3 for details**

## 13. Calculate Full-Scale Concentration

- Transform  $C^*$  to full-scale conditions
- Account for full-scale flow rate and wind speed
- **See Equation 5 for details**

## 14. Optional: Transform to Full Scale or Non-Dimensional

- `full_scale='ms'`: Keep model scale values
- `full_scale='fs'`: Convert all quantities to full scale
- `full_scale='nd'`: Convert all quantities to non-dimensional

## 15. Save Processed Data

- Save time series with all calculated quantities
- Save averaged values (time-averaged statistics)
- Save statistical measures (mean, percentiles, peak-to-mean ratio)
- Create combined CSV file with all measurement points

## 16. Visualization and Analysis

- Plot time series of concentrations
- Generate statistical plots (histograms, PDFs, CDFs)
- Create comparison plots between different measurement locations
- Analyze concentration fluctuations and intermittency
- Plot power spectral density of concentration fluctuations

### 2.2 Example\_Flow\_Analysis (Flow Field Analysis)

This notebook analyzes velocity measurements from wind tunnel experiments:

- **Vertical/Horizontal Wind Profiles:** Analysis of mean velocity and turbulence intensity
- **Turbulence Statistics:** Calculate turbulence intensities ( $I_u$ ,  $I_v$ ,  $I_w$ ), Reynolds stresses
- **Integral Length Scales:** Calculate  $L_{ux}$  from autocorrelation
- **Spectral Analysis:** Power spectral density of velocity fluctuations
- **Convergence Tests:** Bandwidth analysis for different averaging intervals
- **Boundary Layer Fitting:** Power law exponent  $\alpha$  and roughness length  $z_0$

### 2.3 Example\_PointConc\_Analysis\_200325\_Analysis\_deploy

Streamlined deployment version for Jupyter server environments:

- Automatic package installation from GitHub
- Automatic folder structure creation
- Example data download
- Simplified configuration for students/users

### 3 Physical Equations for Concentration Analysis

#### 3.1 Model Scale Mass Flow Rate

The mass flow rate at model scale is calculated accounting for ambient conditions, controller settings, and gas properties:

$$\dot{Q}_{\text{amb}} = \frac{\text{open\_rate} \cdot 10 \cdot f_{\text{gas}} \cdot \dot{Q}_{\text{controller}}}{100} \cdot \frac{T_K \cdot p_{\text{std}}}{p \cdot T_{\text{std},K}} \quad (1)$$

where:

- $\dot{Q}_{\text{amb}}$  – ambient mass flow rate [l/h]
- open\_rate – valve opening [0-10], multiplied by 10 to get percentage
- $f_{\text{gas}}$  – gas correction factor [-] (e.g., for different tracer gases)
- $\dot{Q}_{\text{controller}}$  – maximum flow rate of mass flow controller [l/h]
- $T_K$  – ambient temperature [K]
- $p$  – ambient pressure [Pa]
- $T_{\text{std},K}$  – standard temperature = 273.15 K
- $p_{\text{std}}$  – standard pressure = 101325 Pa

**Note on Mass Flow Controller:** The mass flow controller unit represents the maximum flow rate capacity of the device (e.g., 150 l/h). This value is calibrated at standard conditions and must be corrected for ambient temperature and pressure. The actual flow rate is then determined by the valve opening (open\_rate) which controls the percentage of this maximum flow that passes through.

#### 3.2 Net Concentration

The net concentration is obtained by subtracting the background concentration (slow FID) from the total measured concentration (fast FID):

$$C_{\text{net}} = \text{fast\_FID} - \text{slow\_FID} \quad (2)$$

where:

- $C_{\text{net}}$  – net tracer concentration [ppmV]
- fast\_FID – total concentration measured by fast flame ionization detector [ppmV]
- slow\_FID – background concentration [ppmV]

#### 3.3 Dimensionless Concentration ( $C^*$ )

The dimensionless concentration is the fundamental scaling parameter for comparing measurements at different scales and conditions:

$$C^* = \frac{C_{\text{net}} \cdot \bar{U}_{\text{ref}} \cdot L_{\text{ref}}^2}{\dot{Q}_{\text{amb}}} \quad (3)$$

In computational form (with unit conversions):

$$C^* = \frac{\left(\frac{C_{\text{net}}}{10^6}\right) \cdot \bar{U}_{\text{ref}} \cdot L_{\text{ref}}^2}{\left(\frac{\dot{Q}_{\text{amb}}}{1000 \cdot 3600}\right)} \quad (4)$$

where:

- $C^*$  – dimensionless concentration [-]
- $C_{\text{net}}$  – net concentration [ppmV], converted to volume fraction by  $10^{-6}$
- $\bar{U}_{\text{ref}}$  – mean reference wind speed [m/s]
- $L_{\text{ref}}$  – reference length (model scale) [m]
- $\dot{Q}_{\text{amb}}$  – ambient flow rate [l/h], converted to [m<sup>3</sup>/s] by  $(1000 \cdot 3600)^{-1}$

**Physical Interpretation:**  $C^*$  represents the normalized concentration field. It is scale-invariant, meaning identical  $C^*$  values at different scales indicate similar dispersion physics (assuming geometric and dynamic similarity).

### 3.4 Full-Scale Concentration

To predict concentrations at full scale, the dimensionless concentration is rescaled using full-scale parameters:

$$C_{\text{fs}} = C^* \cdot \frac{\dot{Q}_{\text{fs,vol}}}{L_{\text{ref,fs}}^2 \cdot U_{\text{ref,fs}}} \cdot 10^6 \quad (5)$$

where the full-scale volumetric flow rate is calculated from the mass flow rate:

$$\dot{Q}_{\text{fs,vol}} = \dot{Q}_{\text{fs,mass}} \cdot \frac{R \cdot T_{\text{fs,K}}}{p_{\text{fs}} \cdot M \cdot 10^{-3}} \quad (6)$$

Combined:

$$C_{\text{fs}} = C^* \cdot \frac{\dot{Q}_{\text{fs,mass}} \cdot R \cdot T_{\text{fs,K}}}{p_{\text{fs}} \cdot M \cdot 10^{-3} \cdot L_{\text{ref,fs}}^2 \cdot U_{\text{ref,fs}}} \cdot 10^6 \quad (7)$$

where:

- $C_{\text{fs}}$  – full-scale concentration [ppmV]
- $\dot{Q}_{\text{fs,mass}}$  – full-scale mass flow rate [kg/s]
- $\dot{Q}_{\text{fs,vol}}$  – full-scale volumetric flow rate [m<sup>3</sup>/s]
- $R$  – universal gas constant = 8.3144621 J/(mol·K)
- $T_{\text{fs,K}}$  – full-scale ambient temperature [K]
- $p_{\text{fs}}$  – full-scale ambient pressure [Pa]
- $M$  – molecular weight of air [g/mol], converted to [kg/mol] by  $10^{-3}$
- $L_{\text{ref,fs}}$  – full-scale reference length [m] = scale  $\cdot L_{\text{ref}}$
- $U_{\text{ref,fs}}$  – full-scale reference wind speed [m/s]

**Physical Interpretation:** This transformation accounts for:

- Geometric scaling ( $L_{\text{ref,fs}}^2$  vs  $L_{\text{ref}}^2$ )
- Kinematic scaling ( $U_{\text{ref,fs}}$  vs  $\bar{U}_{\text{ref}}$ )
- Source strength scaling ( $\dot{Q}_{\text{fs,vol}}$  vs  $\dot{Q}_{\text{amb}}$ )
- Thermodynamic conditions (via ideal gas law)



## 4 Parameter File Structure

The ambient conditions are stored in CSV files with the following structure:

### 4.1 Example CSV Format

```
,file1.ts#0,file2.ts#1,file3.ts#2
x_source,0.0,0.0,0.0
y_source,0.0,0.0,0.0
z_source,0.0,0.0,0.0
x_measure,37.5,50.0,62.5
y_measure,0.0,0.0,0.0
z_measure,10.0,10.0,10.0
pressure,101768,101768,101768
temperature,23.0,23.0,23.0
mass_flow_controller,0.300,0.300,0.300
calibration_curve,None,None,None
calibration_factor,1,1,1
scale,200,200,200
ref_length,0.005,0.005,0.005
gas_name,C12,C12,C12
mol_weight,28.97,28.97,28.97
gas_factor,1.0,1.0,1.0
full_scale_wtref,5.0,5.0,5.0
full_scale_flow_rate,0.5,0.5,0.5
full_scale_temp,20,20,20
full_scale_pressure,101325,101325,101325
config_name,test,test,test
```

## 5 Software Usage

### 5.1 Quick Start for Students (Jupyter Server)

1. Upload `Example_PointConc_Analysis_200325_Analysis_deploy.ipynb` and `deploy_config.py` to Jupyter server
2. Open the notebook and run Cell 1 – the software automatically:
  - Installs the `windtunnel` package
  - Creates folder structure (`Data/`, `Results/`)
  - Downloads example data
3. Run subsequent cells to perform analysis
4. To use your own data: Upload files to `Data/InputData/` and `Data/ParameterFiles/`, then modify Cell 2

### 5.2 Local Installation (Windows)

For users who want to work locally with full repository access:

1. Clone or download the repository
2. Run `windows_deploy/setup_WTConc2.bat` (creates virtual environment)
3. Run `windows_deploy/Start_WTConc2_Notebook.bat` to start Jupyter
4. Open any notebook from the browser interface

### 5.3 Manual Installation (Linux/Mac)

```
git clone https://github.com/se04ber/WTSoftwareUtilitiesShare.git
cd WTSoftwareUtilitiesShare
python -m venv WTConc2
source WTConc2/bin/activate # On Windows: WTConc2\Scripts\activate
pip install -r requirements.txt
jupyter notebook
```

## 6 Conclusion

The `WTSoftwareUtilitiesShare` package provides a comprehensive framework for wind tunnel data analysis. The modular structure allows users to:

- Process raw concentration and velocity measurements
- Apply proper scaling transformations
- Perform statistical analysis of turbulent concentration fields
- Generate publication-quality visualizations
- Scale results between model and full scale

The software implements rigorous physical models accounting for:

- Ambient condition corrections (temperature, pressure)

- Mass flow controller calibrations
- Geometric scaling (model to full scale)
- Kinematic scaling (velocity matching)
- Thermodynamic transformations (ideal gas law)

For questions or issues, contact: `Sabrina.ebert@proton.me`

Table 1: Ambient Conditions Parameter File Structure

Parameter	Unit	Description	Importance
x_source	[mm]	X-coordinate of tracer gas source	Defines source position
y_source	[mm]	Y-coordinate of tracer gas source	Defines source position
z_source	[mm]	Z-coordinate of tracer gas source (height)	Defines source position
x_measure	[mm]	X-coordinate of measurement point	Defines measurement location
y_measure	[mm]	Y-coordinate of measurement point	Defines measurement location
z_measure	[mm]	Z-coordinate of measurement point (height)	Defines measurement location
pressure	[Pa]	Ambient atmospheric pressure	Critical for mass flow rate calculation
temperature	[°C]	Ambient temperature	Critical for mass flow rate calculation
mass_flow_controller	[l/h]	Maximum flow rate capacity of controller	Defines maximum source strength
calibration_curve	[-]	Linear calibration coefficient	Corrects controller output
calibration_factor	[-]	Calibration offset	Corrects controller output
scale	[-]	Model to full-scale ratio (e.g., 1:200)	Required for scaling
scaling_factor	[-]	Wind speed correction factor	Adjusts wtref to reference position
ref_length	[m]	Reference length in model scale	Fundamental scaling parameter
ref_height	[m]	Reference height (if applicable)	Boundary layer scaling
gas_name	-	Tracer gas identifier (e.g., C12, N <sub>2</sub> )	Documentation
mol_weight	[g/mol]	Molecular weight of air	Used in ideal gas law
gas_factor	[-]	Correction factor for different gases	Adjusts for tracer properties
full_scale_wtref	[m/s]	Full-scale reference wind speed	Target full-scale velocity
full_scale_flow_rate	[kg/s]	Full-scale source emission rate	Target full-scale source strength
full_scale_temp	[°C]	Full-scale ambient temperature	Full-scale ambient conditions
full_scale_pressure	[Pa]	Full-scale ambient pressure	Full-scale ambient conditions
config_name	-	Configuration identifier	Labels different test conditions